

Bench Protocol – PICO AMC Kit

Additional Resources

- Manuals and Calculator: www.actome.de/resources/downloads
 - PI-Quant Software: pico-bioscience.shinyapps.io/piquant/
-

Notes Before Starting

- Centrifuge: 400 rcf, 5 min, RT (unless otherwise specified)
 - Antibody concentration: – 500 pM for absolute quantification
– 40 pM for qualitative measurement
 - Pipette slowly → avoid air bubbles
 - Recommended: 3 technical replicates per sample
 - For cell analysis: 1 million cells (1×10^4 cells/ μ l)
 - High protein abundance: 6-fold dilution series
 - Low protein abundance: 2-fold dilution series
 - Example protocol describes cultured cells
-

A. Sample Preparation

1. Prepare Buffers & Solutions

Reagent	Preparation Instructions
Additive C (5×)	500 μ l PBS
BSA (5×)	400 μ l PBS
Protease Inhibitor Cocktail (25×)	Dissolve 1 tablet in 2 ml PBS
PIC-PBS	Dissolve 1 tablet in 50 ml PBS
Lysis Buffer (2×)	200 μ l Additive T; 400 μ l Additive C; 200 μ l Additive L; 120 μ l PBS; 80 μ l PIC
Control Buffer (CB)	250 μ l LB-Stock; 100 μ l BSA; 150 μ l PBS
Storage Buffer	20 μ l 10×; 180 μ l PBS (for long-term storage + 8.3 μ l PIC)

2. Cell Lysis


- Harvest cells, wash with ~5 ml PIC-PBS, centrifuge, and discard supernatant.
 - Resuspend in 1 ml PIC-PBS, transfer to a 1.5 ml tube, centrifuge.
 - Count cells, and transfer 1×10^6 cells to a new tube, centrifuge.
 - Resuspend in 100 μ l LB, vortex for 10 s. Lyse for 3 h at 4°C.
 - Sonicate for 5 min (ultrasonic bath, RT).
 - Apply lysate to QIAshredder, centrifuge at 20,000 rcf for 2 min. Transfer flow-through to a new tube.
-

3. Binding Reaction

- Calculate ABX-Mix (Antibodies + LB-Stock + PBS) using the PICO Calculator.

Set up Binding Reaction in a 96-well PCR plate:

NTC	4 μ l LB
ABC	2 μ l CB + 2 μ l ABX (3 replicates)
Sample	2 μ l Lysate + 2 μ l ABX (4 replicates)



- Seal plate \rightarrow vortex for convenient mixture, spin at 1,000 rcf for 30 s \rightarrow incubate overnight at 4°C.

B. dPCR-Mix & Pre-Dilution

- Prepare Master Mix (24-well or 96-well)

Reagents	24-well plate	96-well plate
Ultrapure water	606 μ l	834 μ l
QIAcuity Probe Master Mix	284 μ l	390 μ l
PICO Probe (P8/BL/N6/O7)	45 μ l	62 μ l
PICO Probe (P8/BL/N6/O7)	45 μ l	62 μ l
PICO Probe (P8/BL/N6/O7)	45 μ l	62 μ l
PICO Probe (P8/BL/N6/O7)	45 μ l	62 μ l
Coupling dPCR Mix	36 μ l	50 μ l

Note: If fewer than four PICO Probes are used, fill the missing volume with water.


Pre-Dilution

- Remove adhesive foil from Binding Reaction. Add 36 μ l PBS to each well. Mix well.

For a 24-well PICO Assay:


- Prepare a new 96-well plate for the dilution steps.
- Add PBS to columns 1-3 (for Dilution Series 1, DS1) and columns 5-7 (for DS2). For a 40 pM ABX concentration, add 33 μ l of PBS; for an ABX concentration of 500 pM, add 120 μ l PBS.
- The carry-over volume is always 2 μ l.
- Add 41 μ l of Master Mix to columns 10-12.

e. From DS2, carry over 1 μ l to columns 10-12.



For a 96-well PICO Assay:

- a. Prepare three 96-well plates: Dilution Plate #1, Dilution Plate #2 filled with 39 μ l PBS (for ABX concentration of 40pM) or 143 μ l PBS (for a 500 pM ABX concentration), and Master Mix Plate #3 with 13 μ l Master Mix.
- b. Dilution Plates (#1 and #2): The carry-over volume is 2 μ l.
- c. From Plate #2, carry over 1 μ l to the Master Mix.



- Seal plate according to QIAcuity Manual and start PCR:

Priming: Standard Qiagen Priming Profile

c. PCR:

- 95°C for 2 min
- 40 \times (95°C for 15 s + 58°C for 30 s)

Imaging:

PICO Probe	QIAcuity channel	Integration time	Gain
P8 Probe	FAM, green channel	500 ms	6
BL Probe	HEX, yellow channel	400 ms	6
N6 Probe	TAMRA, orange channel	400 ms	6
O7 Probe	ROX, red channel	300 ms	4

- For data analysis, use the PIQuant Software and User Manual www.actome.de/resources/downloads